



siiora

Preserving Nature's Most Wonderful Creation



Cannulated Screws



www.siiora.com

EDITION - 3

Herbert Screw - Cannulated

2.7mm/ 3.5mm Herbert Screw - Cannulated
- Titanium

Cat No.	Length (mm)
Titanium	
Ti.133.214	14
Ti.133.216	16
Ti.133.218	18
Ti.133.220	20
Ti.133.222	22
Ti.133.224	24
Ti.133.226	26
Ti.133.228	28
Ti.133.230	30
Ti.133.235	35
Ti.133.240	40
Ti.133.245	45
Ti.133.250	50



3.5mm/ 4.5mm Herbert Screw - Cannulated
- Titanium

Cat No.	Length (mm)
Titanium	
Ti.135.314	14
Ti.135.316	16
Ti.135.318	18
Ti.135.320	20
Ti.135.322	22
Ti.135.324	24
Ti.135.326	26
Ti.135.328	28
Ti.135.330	30
Ti.135.335	35
Ti.135.340	40
Ti.135.345	45
Ti.135.350	50



Instruments for Herbert Screw

802.114 Direct Measuring Device for 2.7/3.5mm/
4.5mm Herbert Screws



704.030 Cannulated Screw Driver 2.0mm Tip



704.012 Cannulated Screw Driver 2.5mm Tip



802.009 Cannulated drill Bit Plain Shank - Dia.
2.2mm X 100 mm Length
(for 2.7/3.5mm Herbert Screw)



802.014 Cannulated drill Bit Plain Shank - Dia.
2.5mm X 125mm Length
(for 3.5/4.5mm Herbert Screw)



802.121 Threaded Guide Wire 0.9mm X 225mm Length



802.102 Guide Wire Sleeve for 0.9mm Guide Wire -
Outer Dia. 2.2mm



802.116 Guide Wire Sleeve for 0.9mm Guide Wire -
Outer Dia. 2.5mm



Instruments for Herbert Screw

802.104 Drill Sleeve for 2.2mm Drill Bit



802.106 Drill Sleeve for 2.5mm Drill Bit



802.107 Handle for 0.9 Guide Wire ,2.2mm & 2.5mm Drill Bit Sleeve



400.132 Graphic Aluminium Box of Implants & Instruments For 2.7/3.5 & 3.5/4.5mm Herbert Screws



Instruments for Herbert Cannulated Screw

IMPLANT AND INSTRUMENT SET FOR
2.7/3.5 & 3.5/4.5MM HERBERT SCREW



Cat No. IS.133	Description	Units
TI.133.214	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 14mm	2
TI.133.216	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 16mm	2
TI.133.218	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 18mm	2
TI.133.220	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 20mm	2
TI.133.222	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 22mm	2
TI.133.224	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 24mm	2
TI.133.226	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 26mm	2
TI.133.228	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 28mm	2
TI.133.230	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 30mm	2
TI.133.235	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 35mm	2
TI.133.240	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 40mm	2
TI.133.245	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 45mm	2
TI.133.250	2.7mm/3.5mm Herbert Screw - Cannulated - Titanium Length 50mm	2
TI.135.314	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 14mm	2
TI.135.316	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 16mm	2
TI.135.318	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 18mm	2
TI.135.320	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 20mm	2
TI.135.322	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 22mm	2
TI.135.324	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 24mm	2
TI.135.326	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 26mm	2
TI.135.328	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 28mm	2
TI.135.330	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 30mm	2
TI.135.335	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 35mm	2
TI.135.340	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 40mm	2
TI.135.345	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 45mm	2
TI.135.350	3.5mm/4.5mm Herbert Screw - Cannulated - Titanium Length 50mm	2

Continues.....

Instruments for Herbert Cannulated Screw

Cat No. IS.133	Description	Units
802.114	Direct Measuring Device for 2.7mm/3.5mm/4.5mm Herbert Screws	1
704.030	Cannulated Screw Driver 2.0mm Tip	1
704.012	Cannulated Screw Driver 2.5mm Tip	1
802.009	Cannulated drill Bit Plain Shank - Dia. 2.2mm X 100 mm Length (for 2.7/3.5mm Herbert Screw)	1
802.014	Cannulated drill Bit Plain Shank - Dia. 2.5mm X 125mm Length (for 3.5/4.5mm Herbert Screw)	1
802.121	Threaded Guide Wire 0.9mm X 150mm Length	5
802.102	Guide Wire Sleeve for 0.9mm Guide Wire - Outer Dia. 2.2mm	1
802.116	Guide Wire Sleeve for 0.9mm Guide Wire - Outer Dia. 2.5mm	1
802.104	Drill Sleeve for 2.2mm Drill Bit	1
802.106	Drill Sleeve for 2.5mm Drill Bit	1
802.107	Handle for 0.9 Guide Wire, 2.2mm & 2.5mm Drill Bit Sleeves	1
400.132	Graphic Aluminium Box of Implants & Instruments For 2.7/3.5 & 3.5/4.5mm Herbert Screws	1

Herbert Screw - Cannulated

4.5mm/ 5.5mm Herbert Screw - Cannulated
- Titanium

Cat No.	Length (mm)
Titanium	
Ti.134.230	30
Ti.134.235	35
Ti.134.240	40
Ti.134.245	45
Ti.134.250	50
Ti.134.255	55
Ti.134.260	60
Ti.134.265	65
Ti.134.270	70
Ti.134.275	75
Ti.134.280	80



5.5mm/ 6.5mm Herbert Screw - Cannulated
- Titanium

Cat No.	Length (mm)
Titanium	
TI.136.250	50
TI.136.255	55
TI.136.260	60
TI.136.265	65
TI.136.270	70
TI.136.275	75
TI.136.280	80
TI.136.285	85
TI.136.290	90
TI.136.295	95
TI.136.300	100



Instruments for Herbert Screw

802.115 Direct Measuring Device for 4.5/5.5mm/
6.5mm Herbert Screws



703.016 Cannulated Screw Driver 3.5mm Tip



704.012 Cannulated Screw Driver 2.5mm Tip



802.011 Cannulated Drill Bit Plain Shank - Dia.
3.5mm X 200mm Length
(for 4.5/5.5mm Herbert Screw)



802.012 Cannulated Drill Bit Plain Shank - Dia.
4.0mm X 200mm Length
(for 5.5/6.5mm Herbert Screw)



802.013 Guide Wire 1.5mm X 200mm Length



802.122 Threaded Guide Wire 1.2mm, 200mm Length



802.108 Guide Wire Sleeve for 1.2mm & 1.5mm
Guide Wire - Outer Dia. 3.5mm



Instruments for Herbert Screw

802.117 Guide Wire Sleeve for 1.2mm & 1.5mm
Guide Wire - Outer Dia. 4.0mm



802.110 Drill Sleeve for 3.5mm Drill Bit



802.112 Drill Sleeve for 4.0mm Drill Bit



802.113 Handle for 1.2mm & 1.5mm Guide Wire,
3.5mm & 4.0mm Drill Bit Sleeve



400.133 Graphic Aluminium Box of Implants &
Instruments For 4.5/5.5mm & 5.5/6.5mm
Herbert Screws



Instruments for Herbert Cannulated Screw

IMPLANT AND INSTRUMENT SET FOR
4.5/5.5 & 5.5/6.5MM HERBERT SCREW



Cat No. IS.134	Description	Units
TI.134.230	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 30mm	2
TI.134.235	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 35mm	2
TI.134.240	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 40mm	2
TI.134.245	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 45mm	2
TI.134.250	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 50mm	2
TI.134.255	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 55mm	2
TI.134.260	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 60mm	2
TI.134.265	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 65mm	2
TI.134.270	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 70mm	2
TI.134.275	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 75mm	2
TI.134.280	4.5mm/5.5mm Herbert Screw - Cannulated - Titanium Length 80mm	2
TI.136.250	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 50mm	2
TI.136.255	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 55mm	2
TI.136.260	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 60mm	2
TI.136.265	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 65mm	2
TI.136.270	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 70mm	2
TI.136.275	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 75mm	2
TI.136.280	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 80mm	2
TI.136.285	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 85mm	2
TI.136.290	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 90mm	2
TI.136.295	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 95mm	2
TI.136.300	5.5mm/6.5mm Herbert Screw - Cannulated - Titanium Length 100mm	2

Continues.....

Instruments for Herbert Cannulated Screw

Cat No. IS.134	Description	Units
802.115	Direct Measuring Device for 4.5mm/5.5mm/6.5mm Herbert Screws	1
704.012	Cannulated Screw Driver 2.5mm Tip	1
703.016	Cannulated Screw Driver 3.5mm Tip	1
802.011	Cannulated Drill Bit Plain Shank - Dia. 3.5mm X 200mm Length (for 4.5/5.5mm Herbert Screw)	1
802.012	Cannulated Drill Bit Plain Shank - Dia. 4.0mm X 200mm Length (for 5.5/6.5mm Herbert Screw)	1
802.013	Guide Wire 1.5mm X 200mm Length	5
704.005	Threaded Guide Wire 1.2mm & 150mm	1
802.108	Guide Wire Sleeve for 1.2mm & 1.5mm Guide Wire - Outer Dia. 3.5mm	1
802.117	Guide Wire Sleeve for 1.2mm & 1.5mm Guide Wire - Outer Dia. 4.0mm	1
802.110	Drill Sleeve for 3.5mm Drill Bit	1
802.112	Drill Sleeve for 4.0mm Drill Bit	1
802.113	Handle for 1.2mm & 1.5mm Guide Wire, 3.5mm & 4.0mm Drill Bit Sleeves	1
400.133	Graphic Aluminium Box of Implants & Instruments For 4.5/5.5 & 5.5/6.5mm Herbert Screws	1

Small Cannulated Cancellous Screws

Small Cannulated Cancellous Screw 4.0mm
Short Thread - Stainless Steel & Titanium



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
116.010	10	Ti.116.010
116.012	12	Ti.116.012
116.014	14	Ti.116.014
116.016	16	Ti.116.016
116.018	18	Ti.116.018
116.020	20	Ti.116.020
116.022	22	Ti.116.022
116.024	24	Ti.116.024
116.026	26	Ti.116.026
116.028	28	Ti.116.028
116.030	30	Ti.116.030
116.032	32	Ti.116.032
116.034	34	Ti.116.034
116.036	36	Ti.116.036
116.038	38	Ti.116.038
116.040	40	Ti.116.040
116.042	42	Ti.116.042
116.044	44	Ti.116.044
116.046	46	Ti.116.046
116.048	48	Ti.116.048
116.050	50	Ti.116.050
116.055	55	Ti.116.055
116.060	60	Ti.116.060

Small Cannulated Cancellous Screws

Small Cannulated Cancellous Screw 4.0mm
Half Thread - Stainless Steel & Titanium



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
116.114	14	Ti.116.114
116.116	16	Ti.116.116
116.118	18	Ti.116.118
116.120	20	Ti.116.120
116.122	22	Ti.116.122
116.124	24	Ti.116.124
116.126	26	Ti.116.126
116.128	28	Ti.116.128
116.130	30	Ti.116.130
116.132	32	Ti.116.132
116.134	34	Ti.116.134
116.136	36	Ti.116.136
116.138	38	Ti.116.138
116.140	40	Ti.116.140
116.142	42	Ti.116.142
116.144	44	Ti.116.144
116.146	46	Ti.116.146
116.148	48	Ti.116.148
116.150	50	Ti.116.150
116.152	52	Ti.116.152
116.154	54	Ti.116.154
116.156	56	Ti.116.156
116.158	58	Ti.116.158
116.160	60	Ti.116.160
116.162	62	Ti.116.162
116.164	64	Ti.116.164
116.166	66	Ti.116.166
116.168	68	Ti.116.168

Small Cannulated Cancellous Screws

Small Cannulated Cancellous Screw 4.0mm
Fully Thread - Stainless Steel & Titanium



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
116.210	10	Ti.116.210
116.212	12	Ti.116.212
116.214	14	Ti.116.214
116.216	16	Ti.116.216
116.218	18	Ti.116.218
116.220	20	Ti.116.220
116.222	22	Ti.116.222
116.224	24	Ti.116.224
116.226	26	Ti.116.226
116.228	28	Ti.116.228
116.230	30	Ti.116.230
116.232	32	Ti.116.232
116.234	34	Ti.116.234
116.236	36	Ti.116.236
116.238	38	Ti.116.238
116.240	40	Ti.116.240
116.242	42	Ti.116.242
116.244	44	Ti.116.244
116.246	46	Ti.116.246
116.248	48	Ti.116.248
116.250	50	Ti.116.250
116.255	55	Ti.116.255
116.260	60	Ti.116.260

115.010 Washer for 4.0mm Cancellous Screws
Stainless Steel



Ti.115.010 Washer for 4.0mm Cancellous Screws
Titanium



Instruments for Small Cannulated Cancellous Screw

704.005 Threaded Guide Wire 1.25mm, 150mm



704.026 Cannulated Drill Bit 2.7 / 1.35mm, length 130mm Q.C. End



704.027 Cannulated Drill Bit 3.5/1.35mm Length 130mm Q.C. End



704.008 Cannulated Counter Sink Q.C. 6.0mm



704.009 Cannulated Tap 4.0mm Q.C.



704.010 Double Drill Sleeve 2.7 / 1.25mm



704.011 Drill Guide with Stop 3.5/2.7mm with Two Inserts



704.012 Cannulated Hexagonal Screw Driver - 2.5mm Tip



Instruments for Small Cannulated Cancellous Screw

704.013 Direct Measuring Device - S.S.



704.014 Stylet 1.25mm



704.015 Holding Clip





704.016 Cannulated Quick Coupling Tap Handle with Fiber Handle



307.060 Drill & Tap Sleeve Combined 2.5 /3.5mm



400.088 Graphic Aluminium Instruments Box with One Tray for 4.0mm Cannulated Cancellous Screw



Instruments for Small Cannulated Cancellous Screw

Instruments Set 4.0mm Cannulated
Cancellous Screws

Cat No. IS.022	Description	Units
704.005	Threaded Guide Wire 1.25mm, 150mm	10
704.008	Cannulated Counter Sink Q.C. Dia 6.0mm	1
704.009	Cannulated Tap 4.0mm Q.C.	1
704.010	Double Drill Sleeve 2.7/1.25mm	1
704.011	Drill Guide with Stop 3.5mm/2.7mm with two Inserts	1
704.012	Cannulated Hexagonal Screw Driver - 2.5mm Tip	1
704.013	Direct Measuring Device	1
704.014	Stylet 1.25mm	1
704.015	Holding Clip	1
704.016	Cannulated Quick Coupling Tap Handle with Fibre Handle	1
704.026	Cannulated Drill Bit 2.7/1.35mm, Length 130mm Q.C. End	1
704.027	Cannulated Drill Bit 3.5/1.35mm, Length 130mm Q.C. End	1
307.060	Drill & Tap Sleeve Combined 2.5/3.5mm	1
400.088	Graphic Aluminium Box with One Tray for 4.0mm Cannulated Cancellous Screw	1

400.090 Aluminium Screw Box for 4.0mm Cannulated
Cancellous Screw with Screw Tray & Screw
Holding Forceps (without Implant)



Large Cannulated Cancellous Screw 7.0mm

Large Cannulated Cancellous Screw, 7.0mm
Dia 16mm Thread - Stainless Steel & Titanium



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
111.030	30	Ti.111.030
111.035	35	Ti.111.035
111.040	40	Ti.111.040
111.045	45	Ti.111.045
111.050	50	Ti.111.050
111.055	55	Ti.111.055
111.060	60	Ti.111.060
111.065	65	Ti.111.065
111.070	70	Ti.111.070
111.075	75	Ti.111.075
111.080	80	Ti.111.080
111.085	85	Ti.111.085
111.090	90	Ti.111.090
111.095	95	Ti.111.095
111.100	100	Ti.111.100
111.105	105	Ti.111.105
111.110	110	Ti.111.110
111.115	115	Ti.111.115
111.120	120	Ti.111.120

Large Cannulated Cancellous Screw 7.0mm

Large Cannulated Cancellous Screw, 7.0mm
Dia 32mm Thread - Stainless Steel & Titanium



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
112.030	30	Ti.112.030
112.035	35	Ti.112.035
112.040	40	Ti.112.040
112.045	45	Ti.112.045
112.050	50	Ti.112.050
112.055	55	Ti.112.055
112.060	60	Ti.112.060
112.065	65	Ti.112.065
112.070	70	Ti.112.070
112.075	75	Ti.112.075
112.080	80	Ti.112.080
112.085	85	Ti.112.085
112.090	90	Ti.112.090
112.095	95	Ti.112.095
112.100	100	Ti.112.100
112.105	105	Ti.112.105
112.110	110	Ti.112.110
112.115	115	Ti.112.115
112.120	120	Ti.112.120

Large Cannulated Cancellous Screw 6.5mm

Large Cannulated Cancellous Screw, 6.5mm
Dia 16mm Thread - Stainless Steel



Cat No.	Length (mm)	Cat No.
Stainless Steel		
108.330	30	Ti.108.330
108.335	35	Ti.108.335
108.340	40	Ti.108.340
108.345	45	Ti.108.345
108.350	50	Ti.108.350
108.355	55	Ti.108.355
108.360	60	Ti.108.360
108.365	65	Ti.108.365
108.370	70	Ti.108.370
108.375	75	Ti.108.375
108.380	80	Ti.108.380
108.385	85	Ti.108.385
108.390	90	Ti.108.390
108.395	95	Ti.108.395
108.400	100	Ti.108.400

Large Cannulated Cancellous Screw, 6.5mm
Dia 32mm Thread - Stainless Steel



Cat No.	Length (mm)	Cat No.
Stainless Steel		Titanium
109.430	30	Ti.109.430
109.435	35	Ti.109.435
109.440	40	Ti.109.440
109.445	45	Ti.109.445
109.450	50	Ti.109.450
109.455	55	Ti.109.455
109.460	60	Ti.109.460
109.465	65	Ti.109.465
109.470	70	Ti.109.470
109.475	75	Ti.109.475
109.480	80	Ti.109.480
109.485	85	Ti.109.485
109.490	90	Ti.109.490
109.495	95	Ti.109.495
109.500	100	Ti.109.500

Washer S.S. & Titanium
 For 6.5mm Cancellous &
 For 7.0mm Cannulated Screw

Cat. No.	Cat No.
Stainless Steel	Titanium
115.020	Ti.115.020



Instruments for Large Cannulated Cancellous Screw

703.010 Cannulated Drill Bit 4.5mm



703.011 Cannulated Tap with T - Handle 7.0mm



703.012 Protection Sleeve 4.5mm



703.014 Direct Measuring Device - S.S.
(for 2.0mm Guide Wire)



703.013 Parallel Wire Guide



703.015 Threaded Guide Wire Dia. 2.0mm



703.017 Stylet 2.0mm



703.016 Large Cannulated Hexagonal Screw Driver
3.5mm Tip



6.5mm Large Cannulated Cancellous Instruments

703.010 Cannulated Drill Bit 4.5mm



703.018 Cannulated Tap with T - Handle 6.5mm



703.012 Protection Sleeve 4.5mm



703.014 Direct Measuring Device - S.S.
(for 2.0mm Guide Wire)



703.013 Parallel Wire Guide



703.015 Threaded Guide Wire 2.0mm



703.017 Stylet 2.0mm



703.016 Large Cannulated Hexagonal Screw Driver
3.5mm Tip



Instruments Set for Large Cannulated 6.5mm Cancellous Screw



Cat No.	6.5mm Large Cannulated Cancellous Instrument Set	Units
IS.127		
703.010	Cannulated Drill Bit 4.5mm	1
703.018	Cannulated Tap with T - Handle 6.5mm	1
703.012	Protection Sleeve 4.5mm Dia	1
703.013	Parallel Wire Guide	1
703.014	Direct Measuring Device (for 2.0mm Guide Wire)	1
703.015	Threaded Guide Wire 2.0mm Dia	5
703.016	Large cannulated Hexagonal Screw Driver 3.5mm Tip	1
703.017	Stylet 2.0mm	1
400.089	Graphic Aluminium Instrument Box for 7.0mm Cannulated Cancellous Screw	1

Instruments Box for Large Cannulated Cancellous Screw

400.089 Graphic Aluminium Instrument Box for 7.0mm Cannulated Cancellous Screw



Instrument Set for Large Cannulated Cancellous Screw



Cat No.	7.0mm Large Cannulated Cancellous Instrument Set	Units
IS.020		
703.010	Cannulated Drill Bit 4.5mm	1
703.011	Cannulated Tap with T - Handle 7.0mm	1
703.012	Protection Sleeve 4.5mm Dia	1
703.013	Parallel Wire Guide	1
703.014	Direct Measuring Device (for 2.0mm Guide Wire)	1
703.015	Threaded Guide Wire 2.0mm Dia	5
703.016	Large cannulated Hexagonal Screw Driver 3.5mm Tip	1
703.017	Stylet 2.0mm	1
400.089	Graphic Aluminium Instrument Box for 7.0mm Cannulated Cancellous Screw	1

Cat No.	Large Cannulated Cancellous Implant Set	Units
IS.021		
	Large Cannulated Cancellous Screws 7.0mm Dia, 16mm Thread	
111.030	Length 30mm	3
111.035	Length 35mm	3
111.040	Length 40mm	3
111.045	Length 45mm	3
111.050	Length 50mm	3
111.055	Length 55mm	3
111.060	Length 60mm	3
111.065	Length 65mm	3
111.070	Length 70mm	3
111.075	Length 75mm	3
111.080	Length 80mm	3
111.085	Length 85mm	3
111.090	Length 90mm	3
111.095	Length 95mm	3
111.100	Length 100mm	3
111.105	Length 105mm	2

Instruments and Implants for Large Cannulated Cancellous Screw

Cat No. IS.021	Description	Units
111.110	Length 110mm	2
111.115	Length 115mm	2
111.120	Length 120mm	2
7.0mm Cannulated Cancellous Screws, 32mm Thread		
112.030	Length 30mm	3
112.035	Length 35mm	3
112.040	Length 40mm	3
112.045	Length 45mm	3
112.050	Length 50mm	3
112.055	Length 55mm	3
112.060	Length 60mm	3
112.065	Length 65mm	3
112.070	Length 70mm	3
112.075	Length 75mm	3
112.080	Length 80mm	3
112.085	Length 85mm	3
112.090	Length 90mm	3
112.095	Length 95mm	3
112.100	Length 100mm	3
112.105	Length 105mm	2
112.110	Length 110mm	2
112.115	Length 115mm	2
112.120	Length 120mm	2
115.020	Washer for 6.5 Cancellous and 7.0mm Cannulated Screws	10
400.060	Screw Box for 7.0mm Cannulated Cancellous Screw with Screw Tray. Instrument Tray & Screw Holding Forceps	1

Screw Box

400.060 Screw Box for 7.0mm cannulated Cancellous Screw with Screw Tray, Instrument Tray & Screw Holding Forceps





Siora Surgicals Private Limited

Manufacturer of:
Trauma & Spine Implants



Regd. & Corporate Office

WZ-1, Phool Bagh, Rampura
Delhi - 110035, INDIA



Office Land line:

+91 11 28316116 , 17, 18



Manufacturing Unit:

1792, HSIIDC Indl. Estate
RAI District, Sonapat
Haryana- 131029, INDIA



www.siora.com



impex@siora.net



exports@siora.net